


Renegade


MAN UP PUMP



GENERAL INSTRUCTIONS:

Prior to use, trim pubic hair to ensure a good seal around the base of the penis. If preferred, apply a small amount of lubricant to the head and base of the penis.

1. Massage the penis to encourage blood flow (a semi erect penis can be easier to insert into the pump).
2. Insert the flaccid penis at the entrance of the pump or the semi-erect penis into the pump.
3. Position the pump evenly around the base of the penis.

4. Push the  button and the vacuum will start. You will feel pressure increase, and the penis will be sucked into the chamber. Continue until the penis is fully erect. When your penis is fully erect, you should notice an increase in suction.



5. Press the  button and hold against your body.
6. General recommendations for penile health: pump and create erect penis, leave erect in pump for 3-5 minutes, deflate vacuum pump, wait 1 minute - repeat this 3 times per session. 5-7 sessions/week.
7. Do not use for more than 20-30 minutes at one time.
8. To remove the pump, press the  button and this will reduce the suction.



For warranty information, latest products, and news visit us at:
www.nsnovelties.com